ORTHOFIX® GALAXY FIXATION SYSTEM	EN	ORTHOFIX® GALAXY-FIXATIONSSYSTEM	DA
2		57	
ORTHOFIX® GALAXY FIXATION SYSTEM	IT	ORTHOFIX® GALAXY -KIINNITYSJÄRJESTelmä	FI
13		68	
ORTHOFIX® SYSTEME DE FIXATION GALAXY	FR	ORTHOFIX® GALAXY-FIKSERINGSSYSTEM	NO
24		79	
ORTHOFIX® GALAXY-FIXATIONSSYSTEM	DE	ORTHOFIX® GALAXY FIXATIESYSTEEM	NL
35		90	
ORTHOFIX® SISTEMA DE FIJACIÓN GALAXY	ES		
46			

ORTHOFIX Srl

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Telephone +39 045 6719000 - Fax +39 045 6719380

Instruction leaflet - Please read prior to use

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Telephone 0039-0456719000 - Fax 0039-0456719380

DESCRIPTION & INDICATIONS FOR USE

The Orthofix Galaxy Fixation System consists of a series of components used in conjunction with Orthofix bone gripping elements (bone screws, wires...). External fixation systems are modular, therefore different frame configurations are possible.

The Orthofix Galaxy Fixation System components are not intended to replace normal healthy bone or to withstand the stresses of full weightbearing, particularly in unstable fractures or in the presence of non union, delayed union or incomplete healing. The use of external supports (e.g. walking aids) is recommended as a part of the treatment. The system consists of various modules to be applied in different anatomical sites of the upper and lower limb as well as the pelvis. When used correctly, the Orthofix Galaxy Fixation System maintains limb function, minimizes surgical trauma to anatomical structures, preserves the blood supply and osteogenic potential of the tissues. All Orthofix devices are intended for professional use only. Surgeons who supervise the use of Orthofix devices must have full awareness of orthopaedic fixation procedures as well as adequate understanding of the philosophy of the Orthofix modular system. To promote the proper use of its fixation system, and establish an effective promotional and training tool, Orthofix has developed several manuals or CD-ROM's containing the relevant information (i. e. general philosophy, surgical application, etc.) called "Operative Techniques". These are available in several languages as a complimentary service for surgeons who have adopted the Orthofix system. If you wish to receive a personal copy, please contact Orthofix or its local authorized representative, with a description of the medical device to be used.

INDICATIONS

The Galaxy Fixation System is intended to be used for bone stabilization in trauma and orthopedic procedures, both on adults and all pediatric subgroups excepts newborns as required.

The indications for use include:

- open or closed fractures of the long bones;
- vertically stable pelvic fractures or as a treatment adjunct for vertically unstable pelvic fractures;
- infected and aseptic non-unions;
- joint pathologies/injuries of upper and lower limb, such as:
 - proximal humeral fractures;
 - intra-articular knee, ankle and wrist fractures;
 - delayed treatment of dislocated and stiff elbows;
 - chronic, persistent elbow joint instability;
 - acute elbow joint instability after complex ligament injuries;
 - unstable elbow fractures;
 - additional elbow stabilization of post-operative unstable internal fixation.

The Orthofix Galaxy Wrist external fixator is intended for the following indications:

- intra-articular or extra-articular fractures and dislocations of the wrist with or without soft tissue damage
- polytrauma
- carpal dislocations
- unreduced fractures following conservative treatment
- bone-loss or other reconstructive procedures
- infection

NOTE: The Shoulder Fixation System is intended to be used for proximal humeral fractures where two thirds of the metaphysis is intact.

CONTRAINDICATIONS

The Orthofix Galaxy Fixation System is not designed or sold for any use except as indicated.

Use of the system is contraindicated in the following situations:

- Patients with mental or physiological conditions who are unwilling or incapable of following postoperative care instructions
- Patients with severe osteoporosis*
- Patients with severe, poorly controlled diabetes mellitus
- Patients with compromised vascularity
- Patients with previous infections
- Patients with malignancy in the fracture area
- Patients with neuromuscular deficit or any other conditions that could influence the healing process
- Patients who are HIV positive
- Patients with foreign body sensitivity. Where material sensitivity is suspected, tests should be made prior to implant insertion

* As defined by the World Health Organization: "Bone mineral density of 2.5 standard deviations or more below the mean peak bone mass (average of young, healthy adults) in presence of one or more fragility fractures".

WARNINGS & PRECAUTIONS

1. Fracture stabilization has to be done following correct fracture reduction.
2. The clamp must be closed first manually by turning the metal ring clockwise before locking it firmly by tightening the cam with the Universal T-Wrench or the 5 mm Allen Wrench.
3. Large Clamp (93010), Medium Clamp (93110), Small Clamp (93310), Wrist Module (93350), Wire Locking Clamp (93620), Elbow Hinge (93410), Large-Medium Transition Clamp (93030), Small Multiscrew Clamp-Long (93320), Small Multiscrew Clamp-Short (93330), cannot be disassembled.
4. Before fixator application, make sure that the clamps are fully loosened.
5. Frame stability must be checked intra-operatively before the patient leaves the operating theatre.
6. Compression is never recommended in a fresh fracture.
7. Particular care should be taken to avoid that screws enter the joints or damage the growth plates in children.
8. Depending on the clinical and radiological findings, the surgeon will decide on the number of rods and bone screws needed to achieve the appropriate frame stability.
9. Any device implanted into the patient, such as bone screws and threaded wires, and in general any device which is labelled "single use only": **MUST NOT BE RE-USED.**
10. Screw length and thread length should be selected in accordance with bone and soft tissue dimensions. Excessive penetration of the second cortex by any type of screw should be avoided, because of the risk of soft tissue damage.
11. The maximum diameter of the screw thread should not be greater than one third of the bone diameter (e.g. 6-5 or 6-5.6 mm bone screws for bone diameter greater than 20 mm).
12. For pre-drilled bone screws, pre-drilling with appropriate drill bits and drill guides prior to screw insertion is imperative. Matching grooves on screws and drill bits help the surgeon to use the correct drill bit. Blunt drill bits can cause thermal damage to the bone and should always be discarded.
13. Self-drilling screws with a thread diameter of 5.00 mm or above should never be inserted with a power tool, but always by hand or with a hand drill. Self-drilling screws with smaller thread diameters may be inserted with a power drill at low speed.
14. Transfix pins of 6 mm in diameter are self-drilling and may be inserted with a power drill. These pins are used in association with the Fixator for temporary ligamentotaxis of

the ankle and knee. Orthofix Transfix pins are single use devices and should never be re-used. They are connected to the Galaxy Rods with two Large Clamps.

15. For more stable fixation of a fracture with a fixator, we recommend that the nearest bone screw is applied fairly close to the fracture margin (a minimum of 2 cm is recommended) and that these distances are equal on both sides of the fracture.
16. Do not use MR with non "MR"-marked devices.
17. Appropriate Orthofix instrumentation should be used to insert bone screws correctly.
18. All equipment should be carefully examined prior to use to assure proper working condition. If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.
19. The fixator should be applied at a sufficient distance from the skin to allow for post-operative swelling and for cleaning, remembering that the stability of the system depends upon the bone-fixator distance. If the fixator is sited at a distance of more than 4 cm from the bone the surgeon will decide on the number of rods and bone screws needed to achieve the appropriate frame stability.
20. Components may not be interchangeable between all Orthofix Fixation Systems. Consult individual operative technique guides for interchangeable components.
21. Additional equipment may be required for fixation application and removal such as screw cutters and power drill.
22. Screw and frame integrity should be monitored at regular intervals. To avoid causing injury, the end of threaded wires, bone screws and transfix pins that have been cut, should be protected with special covers.
23. Meticulous screw site hygiene is required.
24. All patients must receive instruction on the use and maintenance of their external fixation frame, and on pin site care.
25. Patients should be instructed to report any adverse or unanticipated effects to the treating surgeon.
26. The fracture site gap should be reassessed periodically during healing, and adjustments to the frame made as necessary. Persistent separation of the fracture ends may lead to delay in bone union.
27. Removal of the device: the surgeon should make the final decision whether a fixation device can be removed.
28. Do not use components of the Orthofix Galaxy Fixation System in conjunction with products of other manufacturers, unless otherwise specified, as the combination is not covered by the necessary validation.
29. **Elbow Distractor**

The Elbow Distractor is intended to be used to distract the joint intra-operatively in case of elbow stiffness.

 - The entity of distraction of the Elbow must be verified under image intensification.
 - It is mandatory to expose the Ulnar nerve prior to distraction.
30. **Shoulder Fixation System**

Shoulder Fixation System: The wire tip should be 5-10 mm from the articular surface of the Humeral head.

 - During wire insertion, use the wire guide to avoid damage of the soft tissues and/or joint impingement. After wire insertion, check joint function.
 - Insert Wires in the safe corridors to avoid damage to anatomical structures.
 - Do not drill Wires into the soft tissues but push them through the skin. Use the drill at low speed when inserting the wires into the bone.
 - 2.5 mm threaded wires are used with the blocking wire clamp.
 - The first threaded wire should be inserted in the centre of the humeral head in order to target its apex.
 - Wires are cylindrical and can be backed up if necessary.
 - Use the dedicated Orthofix instrumentation to insert threaded wires.
 - The skin around the threaded wires (wires sites) must be meticulously cleaned.
31. **Galaxy Wrist**
 - To ensure the correct locking of the multi-screw clamp use always 2 screws and ensure that they are of the same diameter.
 - Depending on the fracture patterns additional fixation techniques may be necessary.

POSSIBLE ADVERSE EFFECTS

1. Nerve or vessel damage resulting from insertion of wires and screws.
2. Superficial or deep bone screw tract infection, osteomyelitis, or septic arthritis, including chronic drainage of bone screw sites after device removal.
3. Oedema or swelling; possible compartment syndrome.
4. Joint contracture, subluxation, dislocation or loss of range of motion.
5. Premature bone consolidation during distraction osteogenesis.
6. Failure of bone to regenerate satisfactorily, development of non union or pseudarthrosis.

7. Fracture of regenerate bone or through bone screw holes after device removal.
8. Loosening or breakage of implants.
9. Bony damage due to inappropriate implants selection.
10. Bone deformity or equinus of the foot.
11. Persistence or recurrence of the initial condition requiring treatment.
12. Reoperation to replace a component or entire frame configuration.
13. Abnormal growth plate development in patients who are skeletally immature.
14. Foreign body reaction to implants or frame components.
15. Tissue necrosis secondary to implants insertion.
16. Pressure on the skin caused by external components when clearance is inadequate.
17. Limb length discrepancy.
18. Excessive operative bleeding.
19. Intrinsic risks associated with anesthesia.
20. Intractable pain.
21. Bone sequestration secondary to rapid drilling of bony cortex with heat build-up and bone necrosis.
22. Vascular disorders including thrombophlebitis, pulmonary embolus, wound hematomas, avascular necrosis.

Warning: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.

Important

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure which require further surgical intervention to remove or replace the external fixation device. Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the external fixation devices are important considerations in the successful utilization of Orthofix external fixation devices by the surgeon. Proper patient selection and the patient's ability to comply with physician instructions and follow prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. If a surgical candidate exhibits any contraindications or is predisposed to any contraindications, DO NOT USE Orthofix Galaxy Fixation devices.

Materials

The Orthofix Galaxy Fixation System is comprised of stainless steel, aluminium alloy, titanium alloy and plastic components. Those components which contact the patient are the percutaneous pins (bone screws), threaded wires, drill bits, guides used during screw insertion. These are manufactured from surgical grade stainless steel. Some of the Orthofix external fixation bone screws (pins) are supplied with a thin, plasma sprayed coating of hydroxyapatite (HA) on the threaded portion of the shaft.

STERILE & NON-STERILE Product

Orthofix provides certain external fixation devices STERILE while others are provided NON-STERILE. Please review the product label to determine the sterility of each device.

Sterile

Devices or kits provided STERILE are labeled as such.

Contents of package are STERILE unless package is opened or damaged. Do not use if package is opened or damaged.

Non-Sterile

Unless otherwise noted, Orthofix external fixation components are provided NON-STERILE. Orthofix recommends that all NON-STERILE components be properly cleaned and sterilized following the recommended cleaning and sterilization procedures. Product integrity and performance are assured only if packaging is undamaged.

Cleaning, Sterilization and Maintenance

Cleaning is an essential pre-requisite to ensure effective disinfection or sterilization. All instruments must be cleaned prior to use and after each re-use of the instruments. Cleaning may be done following the validated cleaning process described below and using washing instruments and machines validated by the hospital.

WARNINGS

- Aluminium based instruments are damaged by alkaline (pH>7) detergents and solutions.
- Anodised coating is damaged by detergents with free halogen ions or sodium hydroxide.
- Detergents and disinfectants with fluoride, chloride, bromide, iodide or hydroxyl ions **MUST NOT** be used.
- **PRODUCTS LABELED FOR SINGLE-USE MUST NOT BE REUSED.**

Limitations and restrictions on reprocessing

- Repeated reprocessing has minimal effect on reusable instruments;
- End of life is normally determined by wear and damage due to use.

Point of use

- It is recommended that instruments are reprocessed as soon as is reasonably practicable following use, putting them inside a sink filled with cold water (< 40°C) for at least 10 minutes and removing gross soiling with a soft cloth or a soft brush.
- Don't use a fixating detergent or hot water as this can cause the fixation of residue which may influence the result of the reprocessing process.

Containment and transportation

- Follow hospital protocols when handling contaminated and bio-hazardous materials;
- All used surgical instruments should be regarded as contaminated. Their handling, collection and transportation should be strictly controlled to minimize any possible risks to patients, personnel and any area of the healthcare facility.

Preparation for decontamination

- Disassemble the devices where applicable. See Orthofix operative techniques for further details.

Cleaning: manual

- 1) Soak the single components in the cleaning solution. Orthofix recommends the use of a pH-neutral enzymatic cleaning agent. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature.
- 2) Brush the single components in the cleaning solution with a surface soft brush until all visible soiling is removed. Use a soft brush to remove residuals from lumens
- 3) Ultrasonic clean the single components in a degassed cleaning solution. Orthofix recommends the use of a pH-neutral or slightly alkaline solution. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature.
- 4) Brush the single components in running tap water
- 5) Carefully hand-dry using absorbent, non-shedding cloth or industrial dryer.
- 6) Rinse the components in sterile or deionized water
- 7) Carefully hand-dry using absorbent, non-shedding cloth or industrial dryer.

Cleaning: Automated

When the devices to be cleaned have lumens or present complexity, a preliminary manual cleaning may be required.

- 1) Place all the instruments into washer baskets.
- 2) Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
- 3) Orthofix recommends the use of a pH-neutral or slightly alkaline solution. Please refer to detergent's manufacturer datasheet for concentration of the solution, required time and temperature.

Orthofix recommends that cycle steps are at least as follows:

1. prewash at low temperatures
2. main wash
3. rinse with demineralized water
4. thermal disinfection
5. drying

4) Rinse the components in sterile or deionized water

5) Carefully hand-dry using absorbent, non-shedding cloth or industrial dryer.

Maintenance, inspection and testing

All instruments and product components should be visually inspected for cleanliness.

If the equipment does not appear visually clean, cleaning steps must be repeated until the products do not result visually clean.

All instruments and product components should be visually inspected for any signs of deterioration that may cause failure in use (such as cracks or damage to surfaces) and functions tested before being sterilized (see detailed Operative Technique Manuals and Instructions for use). Particular attention should be given to:

- Cannulated devices (NB: cannulated drill-bits are single-patient use only).
- Cutting edges: Discard blunt or damaged instruments.
- Hinged instruments: check for smooth movement of hinges without excessive "play".
- Locking mechanisms should be checked for action.
- If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.
- When instruments form part of an assembly, check assembly with matching components.
- Lubricate all parts, except for cam, bush and ball-joint coupling with lubrication oil for medical applications whenever required (see detailed Operative Technique Manuals). The cams and bushes in the ball joints of the Elbow Distractor must be replaced after every use.

Packaging:

- Wrap the tray before sterilization with an approved sterilization wrap or insert it into a rigid sterilization container in order to prevent contamination after sterilization.
- Do not or include additional systems or instruments in the sterilization tray. Sterility cannot be guaranteed if the sterilization tray is overloaded

Sterilization:

- Steam sterilization is recommended. Gas plasma and EtO sterilization should be avoided as they were not validated for Orthofix products;
- Use a validated, properly maintained and calibrated steam sterilizer;
- Do not exceed 140°C (284°F);
- Do not stack trays during sterilization
- Sterilize by steam autoclaving, utilizing a fractioned pre-vacuum cycle or gravity cycle according to the table below:

STEAM STERILIZER TYPE	GRAVITY	PRE-VACUUM
Preconditioning pulses	N/A	4 pulses
Exposure Temperature	132° C (270° F)	132° C (270° F)
Minimum Holding Time	15 minutes	4 minutes
Drying Time	30 minutes	30 minutes

Storage

Store the sterilized instrument in a dry, clean, and dust-free environment at room temperature.

Manufacturer contact:

Please contact your local Orthofix sales representative for further details.

Disclaimer: “The instructions provided above have been validated by Orthofix as being a true description of the preparation of a device for first clinical use or for re-use of multiple use devices. It remains the responsibility of the reprocessor to ensure that the reprocessing, as actually performed using equipment, materials and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process. The cleaning, disinfection and sterilization processes should be adequately recorded. Likewise any deviation by the reprocessor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences and should also be appropriately recorded”.

Note: ANY DEVICE WHICH IS LABELLED “SINGLE USE ONLY” MUST NEVER BE REUSED. ORTHOFIX IS ONLY RESPONSIBLE FOR SAFETY AND EFFECTIVENESS FOR THE FIRST PATIENT USE OF SINGLE USE DEVICES. The institution or practitioner bears full responsibility for any subsequent use of these devices. **CAUTION:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

MRI SAFETY INFORMATION

GALAXY WRIST

Resonance Environment. Non-clinical testing has demonstrated that the Galaxy Wrist Components are MR Conditional. It can be safely scanned under the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla
- Maximum spatial magnetic field gradient of 900-Gauss/cm (90mT/cm)
- Maximum MR System reported-whole-body-averaged specific absorption rate (SAR) of <4.0 W/kg (First Level Controlled Operating Mode).
- No local transmit/receive coils can must be used on the device.
- The Galaxy Fixation System must be entirely outside the MR scanner bore. No part of the Galaxy Fixation System must extend into the MR bore. Therefore MR scanning of body parts where the Galaxy Fixation System is located is Contraindicated.

HEATING INFORMATION

Under the scan conditions defined above, the Galaxy Fixation System frames are expected to produce a maximum temperature rise of 1°C after 15 minutes of continuous scanning.

DISPLACEMENT INFORMATION

The system will not present an additional risk or hazard to a patient in the 3-Tesla 1.5 and MR environment with regard to translational attraction or migration and torque.

GALAXY FIXATION SYSTEM

Galaxy System Fixator Components are labeled MR CONDITIONAL  according to the terminology specified in ASTM F2503 Standard Practice for Marking Medical Devices and Other Items in the Magnetic Resonance Environment.

Non-clinical testing has demonstrated that the Galaxy System Fixator Components is MR Conditional according to the terminology specified in ASTM F2503 Standard Practice for Marking Medical Devices and Other Items in the Magnetic Resonance Environment. Non-clinical testing, done according to ASTM F2052-06, F2213-06, F2182-11, F2119-07, demonstrated that a patient with the Galaxy Fixation System can be safely scanned under the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla
- Maximum spatial magnetic field gradient of 900-Gauss/cm (90mT/cm)
- Maximum whole-body-averaged specific absorption rate (SAR) of 4.0 W/kg in the First Level Controlled Mode for 15 minutes of scanning.
- No local transmit/receive coils can must be used on the device.
- The Galaxy Fixation System must be entirely outside the MR scanner bore.

No part of the Galaxy Fixation System must extend into the MR bore.

Therefore MR scanning of body parts where the Galaxy Fixation System is located is Contraindicated.

Note: All components of Galaxy Fixation System frames must be identified as MR Conditional prior to being placed in or near an MR Environment. The Threaded Wires (93100), the Wire Locking Clamps (93620), and the L Rod (936010) and Semi-Circular Rods (939010, 939020, 939030) are not MR Conditional. Any construct/frame that is using Threaded Wires, the Wire Locking Clamps, the L-Rod and Semi-Circular Rods must therefore be considered as MR Unsafe.

DISPLACEMENT INFORMATION

The system will not present an additional risk or hazard to a patient in the 3-Tesla 1.5 and MR environment with regard to translational attraction or migration and torque.

HEATING INFORMATION

Comprehensive electromagnetic computer modeling and experimental testing was performed on the following systems:

1.5-Tesla/64-MHz: Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner

3-Tesla/128-MHz: Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI, Active-shielded, horizontal field scanner

to determine the worst heating in seven configurations of Orthofix Galaxy Fixation System. From these studies, it is concluded that once the entire external fixation frame is visible outside the MRI bore, the maximum heating is less than 1 degree Celsius. In non-clinical testing the worst scenarios produced the following temperature rises during MRI under

the conditions reported above:

	1.5 Tesla System	3.0 Tesla System
Galaxy Fixation System		
Minutes of scanning	15	15
Calorimetry measured values, whole body averaged SAR (W/kg)	2.2 W/Kg	2.5 W/Kg
Highest temperature Rise less than (°C)	1°C	1°C

Please note that temperature changes reported apply to the designed MR systems and characteristics used. If a different MR system is used, temperature changes may vary but are expected to be low enough for safe scanning as long as all Galaxy System Fixator Components are placed **outside** the MR bore.

MR PATIENT SAFETY

MRI in patients with Galaxy Fixation System can only be performed under these parameters. It is not allowed to scan the Galaxy Fixation System directly. Using other parameters, MRI could result in serious injury to the patient. When the Galaxy Fixation System is used in conjunction with other External Fixation Systems please be advised that this combination has not been tested in the MR environment and therefore higher heating and serious injury to the patient may occur. Because higher in vivo heating cannot be excluded, close patient monitoring and communication with the patient during the scan is required. Immediately abort the scan if the patient reports burning sensation or pain.

Galaxy Fixation System can only be guaranteed for MRI when using the following components to build a frame:

(*the following components are listed in non-sterile configuration. Please consider that the same MRI information and performance are applicable to the same components in gamma-sterile configuration if available (code number preceded by 99- (e.g. 99-93030))

RODS*

Code	Description
932100	Rod 100 mm long, 12 mm diameter
932150	Rod 150 mm long, 12 mm diameter
932200	Rod 200 mm long, 12 mm diameter
932250	Rod 250 mm long, 12 mm diameter
932300	Rod 300 mm long, 12 mm diameter
932350	Rod 350 mm long, 12 mm diameter
932400	Rod 400 mm long, 12 mm diameter
939100	Rod 100 mm long, 9 mm diameter
939150	Rod 150 mm long, 9 mm diameter
939200	Rod 200 mm long, 9 mm diameter

CLAMPS*

Code	Description
93010	Large Clamp
93110	Medium Clamp
93310	Small Clamp
93020	Multi-Screw Clamp
93030	Large-Medium Transition Clamp (Sterile)
93120	Medium Multi-Screw Clamp (Sterile)

939250	Rod 250 mm long, 9 mm diameter
939300	Rod 300 mm long, 9 mm diameter
936060	Rod 60 mm long, 6 mm diameter
936080	Rod 80 mm long, 6 mm diameter
936100	Rod 100 mm long, 6 mm diameter
936120	Rod 120 mm long, 6 mm diameter
936140	Rod 140 mm long, 6 mm diameter
936160	Rod 160 mm long, 6 mm diameter
936180	Rod 180 mm long, 6 mm diameter
936200	Rod 200 mm long, 6 mm diameter

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* Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Orthofix representative if you have questions about the availability of Orthofix products in your area.

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EN All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

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Orthofix Srl